

Fujiiryoki is recognized

World's first in machine inheritance!

Massage chair JP-1100

Specifications

- Main unit dimension (Approx) (During Ottoman receiving) Width 760 x length 1,380 x height 1,250mm (During maximum usage of Arm/leg) Width 880x length 2,000 x height 740mm
- Minimum width before assembly (Approx) 680mm
- Reclining angle (Approx): 110-157 degrees
- Mass (Approx): 78kg
- Power: AC100V (50/60Hz)
- Consumption Power: 130W (Heater: 40W) (Standby power: 0.3W)
- Rated time: 30 Minutes



black (BK)



beige (CS)




brown (BR)

[Mounting-Installation]

- Mounting width Provide the space for reclining and install on a horizontal place. More than 35cm for the foot direction and more than 50cm for the head direction are required.
- Installation Provide the space for reclining and install on a horizontal place. More than 35cm for the foot direction and more than 50cm for the head direction are required.

Mounting width approximately 70cm




[Accessories]

Detachable type back pad/pillow can be easily removed and attached respectively.



[Dimension (mm) \*\*:approx.]



⚠ Safety warning Please read the "Operation Manual" before use and use this product properly.

- The following persons are strongly recommended not to use this product. Certain physical abnormalities may occur. (Including inflammations of subcutaneous tissue), etc.
- Persons prohibited from receiving a massage by a doctor, for example those suffering from thrombosis, severe aneurysm, acute varicose veins, various kinds of dermatitis and skin infections.
- The following persons must consult a medical physician before using this product as certain physical abnormalities may occur. ● Persons using a medical electronic device embedded in the body such as a pacemaker susceptible to electromagnetic disturbances ● People with malignant tumor ● People with heart trouble ● Persons who are pregnant or menstruating, persons thought to be pregnant, or those immediately after childbirth ● Persons who suffer from paresis related to peripheral circulatory disturbance due to diabetes etc. ● Persons with body temperature exceeding 38°C (when strong acute inflammation symptoms, malaise, chills, blood pressure fluctuations are observed, when debilitated) ● Persons with wounds in treatment areas ● Persons who suffer from acute illness (doloric) ● Persons who suffer from osteoporosis ● Persons who suffer from spinal fractures ● Persons suffering from inflammation such as sprains and muscle strain ● Persons who suffer from paresis/thesia ● Persons with abnormalities in the spine ● Persons who suffer from anemia ● Persons who need to massage a part of the body treated medically in the past ● Persons who need bed rest ● Persons receiving medical treatment ● Persons who feel an abnormality on the body apart from these.
- Please do not allow unattended children or persons unable to declare their own intentions to use this chair. Please do not allow infants near the chair.

While using

● The remote control stand can be installed either in the left or right side.

FUJII MEDICAL INSTRUMENTS MFG. CO., LTD.

Please consult with the retailer described in the right side for the FAQs related to the details of this catalogue. Or request to the point of contact of Fujiiryoki customer. For reformation purpose, design and specification might be changed without any prior notification. There might be slight color variations between the printed materials and actual image. Kindly acknowledge. This catalogue description is updated as on October 2015.

FUJIIRYOKI

GOOD DESIGN AWARD Received in 2017-18

Fujiiryoki - Best Model  
Massage Chair JP-1100  
Made in Japan

The best\* chair  
That keeps getting better.



CYBER-RELAX

\*Compared to our other products

# Fujiiryoki has brought together the latest technology to bring you the best massage chair.

Fujiiryoku understands that “each person has different troubles”.

Some of us feel a stiffness from the neck to the upper shoulders, while some of us feel pain inside the shoulder blades.

Everyone suffers from a different problem, be it back strain, lower back pain or tired feet, depending on our lifestyle.

Fujiiryoki's Cyber Relax Series “AS-1100” addresses each and everyone's problems by making use of numerous sophisticated functions.

NEW  
Unique feature

Intensive relaxation  
for each part  
**12 types of  
spot-concentration  
techniques**

NEW

Feel comforted by  
soothing stimuli  
Acupressure on the soles



NEW

Envelop the back  
and soles of the feet  
in warmth with our  
Double Heaters

**Power Stretch**  
It stretches and relaxes the  
entire body to make the  
massage even more effective

Largest in the industry

**Deep Kiwami  
Mecha PRO**  
It targets the deep layers  
of muscles to attack  
the stiffness

Largest number of courses  
in the industry 2

**21 types of  
automatic courses**  
To choose from  
according to various  
individual needs

Unique feature

**Triple-mode  
Air-System**  
3 types of air modes that  
recreate the flexibility  
of the hands

**Flat Position**  
Relaxes the whole body





Unique feature **12 types of spot-concentration techniques**

This offers 7 minute mini courses to intensively relax the stiffness and tiredness typically experienced by people these days due to lack of exercise or desk jobs. The chair is equipped with a total of 12 types of massages – 6 part focused massages, 3 body trunk relaxing techniques and 3 part focused stretches.

Part focused massage **6 types** This offers courses to intensively relieve stiffness and tiredness.

-  **Neck focused massage**  
For the tension in the back of the neck or stiffness in the neck.
-  **Shoulder focused massage**  
Massage targeting the shoulders.
-  **Shoulder blade focused massage**  
For tension and fatigue of the shoulder blades.
-  **Waist focused massage**  
For the tiredness in the waist.
-  **Scott focused massage**  
For the often fatigued buttocks area.
-  **Low focused massage**  
For tiredness in the area behind the knees up to the soles of the feet.



**NEW** **Body Trunk Relaxing Techniques** **3 types** Takes care of the trunk to restore the body's balance.

For stiffness around the shoulder blades, essential for maintaining posture

-  **Back-Trunk relaxing**  
It grips the portion from your shoulders to the upper pelvis and intensively relieves the stiffness around the shoulder blades.
-  **Waist-Trunk relaxing**  
Care targeting the line from your waist up to the buttocks. It locks the pelvis with airbags near the waist and relaxes the muscles from the waist to the buttocks.
-  **Core pelvis stretch**  
For stiffness in the muscles supporting the pelvis. It performs twisting motions with the airbags near the pelvis and gives stretch massages such as extending the muscles of the buttocks.



Smallest number of choices in the industry **21 types of automatic courses**

- 30 minutes courses**
- VIP Course** From warm-up to cool-down, this is a relaxation course in which you can enjoy the popular body heat is often fatigued. We can enjoy luxurious comfort where it feels like you are being massaged and thoroughly by a number of people.
- Summer cooler course** Soothing stimulation of air promotes blood circulation. This is a relaxation course which relieves the stress and fatigue of the entire body.
- Air Relax**
- Slow Stretch** This is a stretch-only course in which you can enjoy the popular stretch movements more slowly and thoroughly.

- 16 minute courses**
- Whole body courses** Effective courses when you want to massage your whole body.
- Whole Body Extreme Relaxation**
- Whole Body Refresh**
- Whole Body Relax**
- Whole Body Sit**
- Neck/shoulder/Extremities Relaxation**
- Neck/shoulder/Extremities Refresh**
- Neck/shoulder/Extremities Relax**
- Whole Body Stretch**

- Neck/shoulder/Extremities Relaxation**
- Neck/shoulder/Extremities Refresh**
- Neck/shoulder/Extremities Relax**
- Waist/Extremities Relaxation**
- Waist/Extremities Refresh**
- Waist/Extremities Relax**
- Waist/Spine/Extremities Relaxation**
- Waist/Spine/Extremities Refresh**
- Waist/Spine/Extremities Relax**

**Alimony courses**  
You can customize all the courses to create your original courses.

- Extreme Whole Body Stretch**
- Extreme Whole Body Relax**
- Whole Body Air**

**Part focused stretches** **3 types**  
This offers courses to stretch the body with stretches specialized for each body part.

-  **Neck stretch**  
It stretches the nape of the neck by fixing the neck with linear balls and extending the entire body with footrests.
-  **Back stretch**  
After stretching your entire body with footrests, it lifts up and arches the back.
-  **Waist stretch**  
It extends the legs and pulls them down while at the same time pushing and stretching the waist up.

**Full Color Touch Panel Remote Control**

An advanced design with a feel of luxury combining visibility and operability.



**Normal Mode**  
Color display and simple touch enables you to select your preferred massage with ease.



**Quick Mode**  
Big letters that are easy to understand. Start your massage with one touch.



© Can switch to English display.

## Air Massage

Wraps your body with 31 air bags and helps you comfortably release the fatigue. To further enhance your comfort, we have installed a feature that lets you adjust the strength for each area.

Air massage is based on repeating the cycle of press and release. With the movements of our airbags that reproduce the power and flexibility of a professional hand massage and the feature that allows you select your preferred level of strength for each body part, you will experience a feeling of blissful vibrancy.



### 1) Shoulder massage

Five-level strength adjustment

Messages by wrapping portion from your shoulders to the upper arms from the outside.



### 2) Arm massage

Three-level strength adjustment

Messages by firmly squeezing the portion from your elbows up to the fingertips.



### 3) Waist/pelvis massage

Five-level strength adjustment

Relaxes a wide range of muscles from around the pelvis to the back of the knee.



### 4) Leg massage

Five-level strength adjustment

Repeats the cycle of compression and release to relax your legs from calves, ankles up to the tops of the feet.



The legs are secured with airbags on the left and right and the area around zusanli line is stimulated. Simultaneously, the airbags at the back gradually apply pressure on the back of the calves to relax them.

## Unique feature Triple-mode air system


Offers 3 air modes that create the flexibility of a hand massage.

It comes with "flow rate adjustment function" for controlling the flow rate and speed of the air. You can choose from 3 modes as per your preference - "Hand kneading mode" that has a varied pace of a hand massage, "Pulse mode" that promotes blood circulation by supplying and releasing air little by little repeatedly, and "Normal mode" that massages the body parts intensively.

Focuses on effectiveness

### Hand Kneading Mode


Applies pressure at a varied pace while maintaining comfortable pressure.



Focuses on promotion of blood circulation


### Pulse Mode

Supports promotion of blood circulation by supplying and releasing air little by little repeatedly.



### Normal Mode

Intensively massage body parts by greatly supplying and releasing air repeatedly.




New

## Foot sole shiatsu

Feet are secured between airbags near the ankles and instep and the rollers near the sole for a shiatsu massage.

### Intensive massage for the soles

We use special rollers that match the arch of the feet's soles. Airbags at the ankles and insteps work in tandem, which holds your feet and gives you a massage that feels like shiatsu. And heaters at the toes help to promote blood circulation.



### 3-level strength adjustment

You can choose a strength for the rollers on the soles from 3 levels by adjusting the strength of the airbags that support your ankles up to the insteps. You can also turn off the rollers if you prefer.



### Unique shape of the rollers

Shape of the rollers offers 2 types of protrusions, high protrusion and gentle protrusion, which stretches and extends the feet from the arch to the outer side of the soles and relaxes them.



Stimulates the inside (high protrusion)

Stimulates the outside (gentle protrusion)



**Soles of our feet are often tired as they support our whole body. It is important to relax the soles (plantar fascia).**


The plantar fascia, extending from the heels to the base of the toes, is a muscle that supports the arch of the sole. It absorbs the impact of walking. It is important to massage and relieve tension from the soles of our feet which always bear the burden of our body.





**Flat Position**

High-quality sleeping posture to help you feel a sense of comfort. Surrender your entire body and relax.



**Power Stretch**

A new way to loosen and relax muscles. Various new stretchers to increase the effectiveness of massages.


**3 types of foot stretchers**



Legs are firmly secured with airbags and pulled downwards. Legs are stretched up by sliding.

**Part-focused stretch**

The body is secured with airbags or knead balls and each part is powerfully stretched by pulling and stretching the legs.



**Perfect fitting footrests**

Electric and springs-type slide mechanism fits one another perfectly.

**1) Footrest slide mechanism**

The space between the main body and the leg rest can be adjusted with a remote control to match the length of your legs.



**2) Springs-type slide**

You can push it out with your feet and it's automatically retracted. It can match it to the length of your legs.



Legs are secured with airbags, stretched up and then pulled down.

**Swing stretch function**

Body is secured with airbags and the backrest and leg part move up and down in tandem to thoroughly stretch the entire body.



\* Incorporated in some courses

\* Compared to our product AS-1000

**Double heaters**

Two types of heaters warm up the muscles around the spinal column, which tend to stiffen, and toes, which tend to get cold.



**Back sensor heater**  
Film heaters built into the back line of the backrest provide effective heating to the area around the spinal column.



**Foot sole heaters**  
Gradually warms the toes, which tend to get cold, and relaxes them.

\* If it is possible to use just the heater \* Pre-set temperature: Approx. 43°C (set temperature and the 'feet' (feet temperature varies).

**Beauty along with healing.**

Sophisticated design with a touch of class, will enhance the quality of your interior space. It comes in 2 colors, elegant and graceful! Black and natural beige.



Beige (S)



Black (Bx)

## MECHANICAL MESSAGE

Largest in the industry

# Deep Kiwami Mecha PRO

Our commitment to the industry's leading Mecha is to remain constantly innovative.

We bring to you knead balls with the largest rate of protrusion in the industry\* (about 12.5 cm). Two knead balls move up/down, right/left and front/back and they can be controlled precisely even at that time. The mecha unit, which combines various kneading techniques, attacks the deep layers of stiffness even more deeply and delicately.



Maximum protruding capacity

Approximately **12.5 cm**

### Unique Feature



#### Double Sensing

It is equipped with 2 sensors to optimize the "depth" and "points" of the message.

The chair is equipped with 2 sensing systems – S-shaped line detection system\* to sense the line of the spine and "shoulder position detection system" to sense the position of the shoulders. It detects the body shape in 3D to adjust the kneading perfectly for each and every person.

Automatically detects the shape of the spine  
Optimizes the "depth" of the message

Back



Shoulders

Automatically detects the position of the shoulders  
Optimizes the "points" of the message

Most numerous in the industry

### Adjust strength up to 12 levels

It allows you to choose a suitable strength according to the stiffness you feel that day.

You can finely adjust the wide amount of protrusion of the knead balls and select your preferred strength to counter any kind of stiffness.

\*Automatic course offers 7 levels, manual course has 12 levels.

Up to **12 levels**  
(When operating manually)



Image Diagram


Unique Feature

### Deep Layer Approach Technique

Seven types of new kneading techniques that work on "deep layers of muscles" which create a high degree of stiffness.

The chair offers "deep layer approach techniques" to relieve the chronic stiffness experienced by people these days more deeply and delicately. Fujityoki brings you knead balls, with the highest rate of protrusion in the industry of about 12.5cm, that can be controlled minutely to work thoroughly on the deep layers of your muscles. From muscles behind the neck supporting the head to the deep muscles around the hips that can cause back ache, it intensively massages each part of your body.

Neck Shiatsu	Loosens and relaxes the "semispinalis capitis muscles" behind the neck that support the head.
Levator scapulae muscle relax	Kneads and relaxes the muscles that pull down the shoulder blades. Effective for stiffness in shoulders caused by bad posture.
Rhomboid muscle relax	Increases strength and relaxes the muscles that pull the scapular blades back when you stretch your chest (trapezius muscles).
Spinal muscle stretch	Presses and stretches the muscles around the spinal column supporting the torso.
Multifidus lumborum relax	Relaxes the deep layers of muscles that cause pain in the back or the waist.
Psoas major muscles shiatsu	Thoroughly presses and relaxes the muscles supporting the pelvis, like shiatsu.
Buttocks muscles relax	Massages deep layers of muscles around the buttocks that cause lower back pain.




Top class in the industry

### 41 types, 85 kinds of kneading techniques

Offers a variety of specialized techniques targeting all the areas such as the Neck, Shoulders, Back, Waist, Hips.

Various part focused massages utilizing the kneading balls to their fullest extent. It comes with 85 top class kneading techniques, including spot-specific movements. It offers personalized special techniques such as "loop kneading" and "3D wavelike" techniques.



Kneading balls grasp and stretch the neck which promotes blood circulation around the neck.  
Stimulates the shoulders deeply from the top to relieve the stiffness in the trapezius from fatigue and neuralgia or muscle pain.  
Effectively relaxes the waist and buttocks from fatigue and neuralgia or muscle pain.

Knead up	Knead down	Roll back & forward	Rolling	Wavelike	Palmar knead up
Palmar knead down	Neck relax	Extreme knead	Shoulder blade relax	Shoulder blade press	Loop knead down
Loop tapping	Loop wavelike	Stretch and shiatsu	Shiatsu Shiatsu	Wave rolling	Push and knead
Rolling	Para spinal	Psoas major muscle stretch	Extreme waist tapping	Stretch	Stretch and wavelike
Stretch knead up	3D knead up	3D tapping	3D wavelike	Buttock	Levator scapulae muscle relax
Palmar/wavelike					Rhomboid muscle relax